

For Research Use Only

CoraLite® Plus 750 Mouse IgM Isotype Control (MOPC-104E)



Catalog Number: CL750-65243

Basic Information

Catalog Number:

CL750-65243

Size:

100ug, 200 µg/ml

Source:

Mouse

Isotype:

IgM, lambda

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

N/A

CloneNo.:

MOPC-104E

Excitation/Emission maxima
wavelengths:

755 nm / 780 nm

Applications

Tested Applications:

FC

Species Specificity:

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 human whole blood cells were surface stained with 0.2 μg CoraLite® Plus 750 Anti-Mouse IgM Isotype Control (red), 5 μl CoraLite 750 Plus Anti-Human STRO-1 (green), or unstained (blue). Cells were not fixed. CD45dim, CD34+ lymphocytes were gated.