

For Research Use Only

# CoraLite® Plus 750 Anti-Human IL3RA/CD123 (HI12H7)

Catalog Number: CL750-65256



## Basic Information

### Catalog Number:

CL750-65256

### Size:

100 tests, 5 ul/test

### Source:

Mouse

### Isotype:

IgG1, kappa

### GenBank Accession Number:

BC035407

### GeneID (NCBI):

3563

### ENSEMBL Gene ID:

ENSG00000185291

### Full Name:

interleukin 3 receptor, alpha (low affinity)

### Calculated MW:

378 aa, 43 kDa

### Purification Method:

Affinity purification

### CloneNo.:

HI12H7

### Excitation/Emission maxima wavelengths:

755 nm / 780 nm

## Applications

### Tested Applications:

FC

### Species Specificity:

human

## Background Information

CD123, also referred to as interleukin-3 receptor subunit alpha (IL-3R  $\alpha$ ), belongs to the cytokine receptor family. It is a glycoprotein containing an extracellular domain, involving a predicted Ig-like domain, two FnIII domains, a transmembrane domain and an intracellular domain (PMID: 24513123). CD123 is expressed on hematopoietic progenitor cells, plasmacytoid dendritic cells, basophils, monocytes, some B cells, and in a variety of hematologic malignancies (PMID: 10527461; 34855461). IL-3 initially binds to CD123 and subsequently recruits the beta chain (CD131) to form the high-affinity receptor, resulting in the activation of a series of signaling pathways including the JAK/STAT, Ras-MAPK, and phosphatidylinositol 3-kinase pathways (PMID: 29296749). IL-3R signaling regulates the proliferation, survival, and differentiation of hematopoietic cells.

## Storage

### Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> human PBMCs were surface co-stained with PE Anti-Human CD45RO and 5 ul CoraLite® Plus 750 Anti-Human CD123 (CL750-65256, Clone:HI 12H7), or Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.