

For Research Use Only

# CoraLite® Plus 750 Anti-Mouse GL7 (GL-7)



Catalog Number: **CL750-65261**

## Basic Information

Catalog Number:

CL750-65261

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgM, kappa

GenBank Accession Number:

GeneID (NCBI):

104230

Full Name:

lymphocyte antigen 77

Purification Method:

Affinity purification

CloneNo.:

GL-7

Excitation/Emission maxima wavelengths:

755 nm / 780 nm

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> unstimulated or anti-CD3/CD28 stimulated mouse splenocytes were surface co-stained with CoraLite® Plus 488 Anti-Mouse CD3 (17A2) and 0.5 ug CoraLite® Plus 750 Anti-Mouse GL7 (CL750-65261, Clone:GL-7). Cells were not fixed.