

For Research Use Only

Hexylene glycol

Catalog Number: CM00293

产品信息

Catalog Number:
CM00293

CAS号:
107-41-5


分子式:
 $C_6H_{14}O_2$

主要靶点:
Others

主要通路:
其他

分子量:
118.17

溶解度:
H₂O:100 mg/mL (846.24 mM), Sonification is recommended., DMSO:100 mg/mL (846.24 mM)



描述

Hexylene glycol (Diolane), colorless liquid, is a chiral diol. It is an industrial compound produced from diacetone alcohol by hydrogenation.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years