

Catalog Number: CM00297

产品信息

Catalog Number:
CM00297

CAS号:
482-35-9


分子式:
C₂₁H₂₀O₁₂

主要靶点:
NO Synthase|Wnt/beta-catenin|NF- κ B

主要通路:
NF- κ B信号通路|干细胞|免疫与炎症|细胞骨架

分子量:
464.37

溶解度:
DMSO:83 mg/mL (178.74 mM)



描述

1. Isoquercitrin is a potential stimulator of bone mineralization used for prophylaxis of osteoporotic disorders. 2. Isoquercitrin may be as a potential therapeutic agent against neurodegeneration in Parkinson's disease. 3. Isoquercitrin is an inhibitor of Wnt/ β -catenin and may be as a potential novel anti-tumoral agent, such as against human pancreatic liver cancer related to opioid receptors and to the activation of the mitogen-activated protein kinase (MAPK) signalling pathway..

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years