

For Research Use Only

DSG Crosslinker

Catalog Number: CM02037

产品信息

Catalog Number:
CM02037


分子量:
326.26

CAS号:
79642-50-5

分子式:
 $C_{13}H_{14}N_2O_8$

主要靶点:
Others

主要通路:
其他



体外活性

ADCs consist of an antibody attached to an ADC cytotoxin through an ADC linker.

描述

DSG Crosslinker is a cleavable ADC linker used in antibody-drug conjugates (ADCs) synthesis.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.