

For Research Use Only

Fmoc-NH-PEG4-CH₂CH₂COOH

Catalog Number: CM02058

产品信息

Catalog Number:
CM02058

CAS号:
557756-85-1

分子式:
C₂₆H₃₃NO₈

主要靶点:
ADC Linker

主要通路:
ADC

分子量:
487.5421

MDL NO:
MFCD06656472

Pubchem ID:
2756186



描述

Fmoc-NH-PEG4-CH₂CH₂COOH is a cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs).

储存

储存条件:

粉末	-20°C	3年
----	-------	----

液体	-80°C	1年
----	-------	----

运输条件:

Shipped in cold pack

参考文献

- 1 Cai P, et al. Synthesis of C3-Alkylated Indoles on DNA via Indolyl Alcohol Formation Followed by Metal-Free Transfer Hydrogenation. Org Lett. 2019 Aug 14.