

Catalog Number: CM02245

产品信息

Catalog Number:
CM02245

CAS号:
160743-62-4


分子式:
 $C_{34}H_{69}NO_6$

主要靶点:
Others

主要通路:
其他

分子量:
587.927

溶解度:
DMSO:100 mg/mL (Need
ultrasonic)



描述

DMG-PEG2000 is a lipid excipient that has been used in combination with other lipids in the formation of lipid nanoparticles. Formulations containing DMG-PEG2000 have been used in the development of lipid nanoparticles for the delivery of mRNA-based vaccines.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years