

For Research Use Only

# Solvent Blue 35

Catalog Number: CM02252

## 产品信息

**Catalog Number:**  
CM02252

**CAS号:**  
17354-14-2


**分子式:**  
 $C_{22}H_{26}N_2O_2$

**主要靶点:**  
Others

**主要通路:**  
其他

**分子量:**  
350.45

**溶解度:**  
DMSO:40 mg/mL (114.14 mM), Sonification is recommended.



## 描述

Solvent Blue 35 (Sudan Blue II) is a dye used for coloring hydrocarbon-based and alcoholic solvents. It is used for staining triglycerides in animal tissues.

## 储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years