

For Research Use Only

3-(4-Acetylphenyl)-2-Aminopropanoic Acid HCl



Catalog Number: CM02492

产品信息

Catalog Number:
CM02492

CAS号:
1360436-95-8

分子式:
C₁₁H₁₄ClNO₃


主要靶点:
ADC Linker

主要通路:
ADC

分子量:
243.6868

MDL NO:
MFCD08436130

Pubchem ID:
46738304.0



描述

3-(4-acetylphenyl)-2-aminopropanoic acid is a PEGn linker for antibody-drug-conjugation (ADC).

储存

储存条件:

粉末	-20°C	3年
液体	-80°C	1年

运输条件:

Shipped in cold pack

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.