

For Research Use Only

Fmoc-NH-PEG6-CH₂COOH

Catalog Number: CM02502

产品信息

Catalog Number:
CM02502

CAS号:
437655-96-4


分子式:
C₂₉H₃₉NO₁₀

主要靶点:
ADC Linker

主要通路:
ADC

分子量:
561.6207

MDL NO:
MFCD26142983



描述

Fmoc-NH-PEG6-CH₂COOH is a cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs).

储存

储存条件:

粉末	-20°C	3年
----	-------	----

液体	-80°C	1年
----	-------	----

运输条件:

Shipped in cold pack

参考文献

- 1 Michael A, et al. Synthesis of Bifunctional Integrin - Binding Peptides Containing PEG Spacers of Defined Length for Non - Viral Gene Delivery. Volume2008, Issue17.