


Catalog Number: CM03150

产品信息

Catalog Number:
CM03150CAS号:
1461-15-0分子式:
 $C_{30}H_{26}N_2O_{13}$ 主要靶点:
Others主要通路:
其他分子量:
622.53溶解度:
H₂O:5 mg/mL (8.03
mM),ultrasonic and adjust pH to
12 with NaOH,DMSO:100 mg/mL
(160.63 mM),Need ultrasound

描述

Calcein is a fluorescent dye and self-quenching probe that used as an indicator of lipid vesicle leakage, and a complexometric indicator for titration of calcium ions with EDTA, and for fluorometric determination of calcium.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years