

For Research Use Only

Thalidomide-O-C6-NH₂ hydrochloride



Catalog Number: CM03875

产品信息

Catalog Number:
CM03875

CAS号:
2245697-88-3


分子式:
C₁₉H₂₄ClN₃O₅

主要靶点:
E3 Ligase Ligand-Linker
Conjugate

主要通路:
PROTAC

分子量:
409.86

溶解度:
DMSO:83.33 mg/mL (195.19
mM), Sonification is
recommended., H₂O:25 mg/mL
(61.00 mM), Sonification is
recommended.



描述

Thalidomide-O-C6-NH₂ hydrochloride is a synthesized E3 ligase ligand-linker conjugate.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.