


Catalog Number: CM04315

产品信息

Catalog Number:
CM04315**CAS号:**
36951-72-1**分子式:**
 $C_{72}H_{66}N_8O_{12}S_4$ **主要靶点:**
Telomerase|DNA/RNA Synthesis**主要通路:**
DNA损伤和修复|细胞周期**分子量:**
1363.6**溶解度:**
DMSO:10 mg/mL (7.33 mM), Sonification is recommended., H₂O:6.2 mg/mL (4.55 mM), Sonification is recommended.

描述

TMPyP4 tosylate is a quadruplex-specific ligand and is a telomerase inhibitor with antitumor effects in osteosarcoma cell lines.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years