

For Research Use Only

Givinostat hydrochloride monohydrate



Catalog Number: CM04410

产品信息

Catalog Number:
CM04410

CAS号:
732302-99-7


分子式:
 $C_{24}H_{27}N_3O_4 \cdot HCl \cdot H_2O$

主要靶点:
HDAC

主要通路:
DNA损伤和修复|表观遗传

分子量:
475.97

溶解度:
Ethanol:1 mg/mL(2.1 mM),DMSO:88 mg/mL (184.9 mM),H₂O:<1 mg/mL



靶点活性

HD1-B:7.5 nM|HD1-A:16 nM|HD2:10 nM

描述

Givinostat hydrochloride monohydrate is an HDAC inhibitor.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.