


Catalog Number: CM04749

## 产品信息

Catalog Number:  
CM04749CAS号:  
36703-88-5分子式:  
 $C_{52}H_{78}N_{10}O_{17}$ 主要靶点:  
HIV Protease|IL Receptor主要通路:  
免疫与炎症|微生物学|蛋白酶体分子量:  
1115.23溶解度:  
DMSO:89.7 mM,H<sub>2</sub>O:89.7 mM

## 描述

Inosine pranobex (Delimmun) has a number of immunomodulatory effects, including inducing T-lymphocyte differentiation, augmenting macrophage and NK cell functions, and stimulating IL-2 production. Through these effects and others, isoprinosine shows antiviral activity and has applications against subacute sclerosing panencephalitis. It also restores depressed immune responses in Y patients after irradiation and reduces the incidence of infection in leukemia patients undergoing chemotherapy. It is most commonly used to treat the rare measles complication subacute sclerosing panencephalitis in conjunction with intrathecal interferon therapy.

## 储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years