

Catalog Number: CM06278

产品信息

Catalog Number:
CM06278

CAS号:
153559-46-7


分子式:
C₂₃H₂₆O₃

主要靶点:
Retinoid Receptor|Autophagy

主要通路:
代谢|自噬

分子量:
350.45

溶解度:
DMSO:49 mg/mL (139.82 mM)



靶点活性

RXR α :330 nM (EC50)|RXR γ :260 nM (EC50)|RXR β :200 nM (EC50)

体外活性

LG-100064 shows no effect on RAR α / β / γ (EC50, >10000 nM)[1]. LG-100064 (3-Methyl-TTNCB) is a retinoid-X-receptor (RXR) agonist, with EC50s of 330 nM, 200 nM, and 260 nM for RXR α , RXR β and RXR γ .

储存

Powder: -20°C for 3 years | In solvent: -80°C for 1 year | Shipping with blue ice.