

For Research Use Only


Boc-NH-PEG4-CH₂COOH

Catalog Number: CM11506

产品信息

Catalog Number:
CM11506

分子量:
351.39



CAS号:
876345-13-0

分子式:
C₁₅H₂₉NO₈

主要靶点:
PROTAC Linker

主要通路:
PROTAC

描述

Boc-NH-PEG4-CH₂COOH is a cleavable ADC linker used as a linker for antibody-drug conjugates (ADC)[1].

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years