

For Research Use Only

Poloxamer 407



Catalog Number: CM15882

产品信息

Catalog Number:
CM15882

CAS号:
9003-11-6

分子式:

主要靶点:
Others

主要通路:
其他

溶解度:

H2O:100 mg/mL, Sonification is recommended.

Poloxamer 407

描述

Poloxamer 407 is a nonionic surfactant.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.