

For Research Use Only

Activated Protein C (390-404), human acetate



Catalog Number: CM16077

产品信息

Catalog Number:
CM16077

分子量:
1960.19

CAS No.:

CAS号:

分子式:

$C_{93}H_{134}N_{22}O_{25}$

主要靶点:

APC

主要通路:

细胞周期

储存

keep away from moisture | Powder: -20°C for 3 years | In solvent: -80°C for 1 year | Shipping with blue ice.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.