

For Research Use Only

Adropin (34-76) (human, mouse, rat)



Catalog Number: CM16310

产品信息

Catalog Number:
CM16310

分子量:
4499.82

CAS号:
1802086-30-1

分子式:
 $C_{190}H_{293}N_{55}O_{68}S_2$

主要靶点:

L-Proline, L-cysteinyl-L-histidyl-L-seryl-L-arginyl-L-seryl-L-alanyl-L- α -aspartyl-L-valyl-L- α -aspartyl-L-seryl-L-leucyl-L-seryl-L- α -glutamyl-L-seryl-L-seryl-L-prolyl-L-asparaginyl-L-seryl-L-seryl-L-prolylglycyl-L-prolyl-L-cysteinyl-L-prolyl-L- α -glutamyl-L-lysyl-L-alanyl-L-prolyl-L-prolyl-L-prolyl-L-glutamyl-L-lysyl-L-prolyl-L-seryl-L-histidyl-L- α -glutamylglycyl-L-seryl-L-tyrosyl-L-leucyl-L-leucyl-L-glutamyl-, cyclic (1-23)-disulfide

描述

Adropin (34-76) (human, mouse, rat) is a peptide that attenuates liver fibrosis and insulin resistance independent of obesity or food intake.

储存

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.