

Catalog Number: CM20313

产品信息

Catalog Number:
CM20313

CAS号:
T11855


分子式:
 $C_{205}H_{366}N_3O_{117}P_5$

主要靶点:
Others

主要通路:
其他

分子量:
4899.92

溶解度:
H2O:5 mg/mL (Need ultrasonic)



描述

Lipopolysaccharide is a highly immunogenic antigen which can enhance immune responses. Lipopolysaccharides is an endotoxin derived from the outer leaflet of the outer membrane of Gram-negative bacteria. Lipopolysaccharide consists of an antigen O-specific chain, a core oligosaccharide and lipid A.

储存

Powder: -20°C for 3 years | In solvent: -80°C for 2 years