

For Research Use Only

PerCP Anti-Mouse CD8a Rabbit Recombinant Antibody

Catalog Number: CP-98472-2



Basic Information

Catalog Number: CP-98472-2	GenBank Accession Number: NM_001081110.2	Purification Method: Protein A purification
Concentration: 100ug, 200 ug/ml	GeneID (NCBI): 12525	CloneNo.: 242879C10
Source: Rabbit	UNIPROT ID: P01731	Recommended Dilutions: FC: 0.25 ug per 10 ⁶ cells in a 100 µl suspension
Isotype: IgG	Full Name: CD8 antigen, alpha chain	Excitation/Emission maxima wavelengths: 482 nm / 677 nm
Immunogen Catalog Number: EG3265	Calculated MW: 27 kDa	

Applications

Tested Applications: FC	Positive Controls: FC : mouse splenocytes,
Species Specificity: mouse	

Background Information

CD8 is a transmembrane glycoprotein composed of two disulfide-linked chains. It can be present as a homodimer of CD8 α or as a heterodimer of CD8 α and CD8 β (PMID: 3264320; 8253791). CD8 is found on most thymocytes. The majority of class I-restricted T cells express mostly the CD8 α β heterodimer while CD8 α α homodimers alone have been found on some gut intraepithelial T cells, on some T cell receptor (TCR) γ δ T cells and on NK cells (PMID: 2111591; 1831127; 8420975). CD8 acts as a co-receptor that binds to MHC class-I and participates in cytotoxic T cell activation (PMID: 8499079). During T cell development, CD8 is required for positive selection of CD4-/CD8+ T cells (PMID: 1968084).

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 0.09% sodium azide and 0.5% BSA, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ mouse splenocytes were surface stained with CoraLite® Plus 647 Anti-Mouse CD3 and 0.25 ug PerCP Anti-Mouse CD8a Rabbit RecAb (CP-98472-2, Clone:242879C10), or unstained cells. Cells were not fixed.



1x10⁶ mouse splenocytes were surface stained with 0.25 ug PerCP Anti-Mouse CD8a Rabbit RecAb (CP-98472-2, Clone:242879C10)(red), or unstained cells(blue). Cells were not fixed.