

For Research Use Only

PerCP-Cyanine5.5 Mouse IgG1 Isotype Control (MOPC-21)

Catalog Number: CPY5-65124



Basic Information

Catalog Number:

CPY5-65124

Size:

100ug, 200 ug/ml

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

GeneID (NCBI):

16017

UNIPROT ID:

P01868

Full Name:

immunoglobulin heavy constant gamma 1 (G1m marker)

Purification Method:

Unknow

CloneNo.:

MOPC-21

Excitation/Emission maxima wavelengths:

482 nm / 690nm

Applications

Tested Applications:

FC

Species Specificity:

n/a

Background Information

The MOPC-21 immunoglobulin is useful as an isotype-matched control. The MOPC-21 immunoglobulin has an unknown binding specificity and is used as an isotype control for mouse IgG1 antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human PBMCs were surface stained with 0.1 ug PerCP-Cyanine5.5 Mouse IgG1 Isotype Control (MOPC-21) (CPY5-65124, Clone: MOPC-21) (red), or PerCP-Cyanine5.5 Anti-Human CD45 (HI30) (CPY5-65109, Clone: HI30) (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.