

For Research Use Only.  
Not For Use In Diagnostics.

# PerCP-Cyanine5.5-conjugated Streptavidin Rabbit RecAb

Catalog Number: CPY5-PF00030



---

## Information

Catalog Number:

CPY5-PF00030

Host:

Applications:

FC, IF

Reactivity:

Physical State:

Liquid

Conjugation:

PerCP-Cyanine5.5

---

## Recommended Dilutions

1:200-1:1000 for flow cytometry and immunofluorescence

---

## Storage

Storage:

2-8°C for one year.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA, pH 7.4.

Aliquoting is unnecessary

---

## Purity

---

For technical support and original validation data for this product please contact:

T: 4006900926


E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> human PBMCs were surface stained with Biotin Anti-Human CD8 (UCHT4) Mouse IgG2a Recombinant Antibody (Biotin-65559, Clone: UCHT4), washed with PBS, and stained with PerCP-Cyanine5.5 Streptavidin Recombinant Antibody (CPY5-PF00030) and analyzed by flow cytometry. Cells were not fixed. Lymphocytes were gated.