

For Research Use Only

# Cardinal Red™ Anti-Human CD14 (UCHM-1)



Catalog Number: CR-65056

## Basic Information

<b>Catalog Number:</b> CR-65056	<b>GenBank Accession Number:</b> BC010507	<b>Purification Method:</b> Affinity purification
<b>Size:</b> 100tests, 5 µl/test	<b>GeneID (NCBI):</b> 929	<b>CloneNo.:</b> UCHM-1
<b>Source:</b> Mouse	<b>ENSEMBL Gene ID:</b> ENSG00000170458	<b>Excitation/Emission maxima wavelengths:</b> 592 nm / 611 nm
<b>Isotype:</b> IgG2a	<b>UNIPROT ID:</b> P08571	
	<b>Full Name:</b> CD14 molecule	
	<b>Calculated MW:</b> 375 aa, 40 kDa	

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
Human

## Background Information

CD14 is a 50-55 kDa glycosylphosphatidylinositol-anchored glycoprotein preferentially expressed on monocytes and macrophages, and at lower levels on granulocytes (PMID: 3385210; 2462937; 7685797). CD14 can also exist as a soluble protein. CD14 acts as a co-receptor for bacterial liposaccharides (LPS) (PMID: 1698311). It plays a major role in the inflammatory response of monocytes to LPS.

## Storage

**Storage:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> human PBMCs were surface stained with 5 ul Cardinal Red™ Anti-Human CD14 (CR-65056, Clone:UCHM-1) or unstained. Cells were not fixed. Monocytes were gated.



1X10<sup>6</sup> human PBMCs were surface stained with 5 ul Cardinal Red™ Anti-Human CD14 (CR-65056, Clone:UCHM-1) or unstained. Cells were not fixed. Monocytes were gated.