

For Research Use Only

Cardinal Red™ Mouse IgG2a isotype control (11A1B2)



Catalog Number: CR-66360-2

Basic Information

Catalog Number:

CR-66360-2

Size:

100ug, 0.2 mg/mL

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

11A1B2

Excitation/Emission maxima
wavelengths:

592 nm / 611 nm

Applications

Tested Applications:

FC

Species Specificity:

Megathura crenulata

Background Information

This antibody is used as the isotype control of mouse IgG2a. It can be used as the isotype control in Flow Cytometry, Immuno-Precipitation and other experiments. The immunogen of this antibody is Keyhole Limpet Hemocyanin (KLH).

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 human PBMCs were surface stained with 0.2 μ g Cardinal Red™ Anti-Mouse IgG2a isotype control (CR-66360-2, Clone:11A1B2) (red), Cardinal Red™ Anti-Human CD3 (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.