

Recombinant Mouse IL-1 alpha protein (His Tag)(HPLC verified)

Catalog Number: Eg0263

Basic Information

Species:
Mouse**Purity:**
>90 %, SDS-PAGE
>90 %, SEC-HPLC**Tag:**
His Tag

Technical Specifications

Purity:
>90 %, SDS-PAGE
>90 %, SEC-HPLC**Endotoxin Level:**
<0.1 EU/ µg protein, LAL method**Source:**
HEK293-derived Mouse IL-1 alpha protein Met+Ser115-Ser270 (Accession# P01582) with a His tag at the C-terminus.**GeneID:**
16175**Accession:**
P01582**Predicted Molecular Mass:**
21.9 kDa**SDS-PAGE:**
24-32 kDa, reducing (R) conditions**Formulation:**
Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Not tested

Storage and Shipping

Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

Interleukin-1 α (IL-1 α) is a protein that plays a crucial role in the regulation of a number of key cellular processes. IL-1 α is a cytokine and a potent mediator of body's response to inflammation, microbial invasion, tissue injury and immunological response. In addition, recent studies suggest that IL-1 α also has a role in wound healing, rheumatoid arthritis, Alzheimer's disease and tumor growth. IL-1 α functions as an alarmin during the process of sterile inflammation. It serves as a signal from dying cells to potentiate the inflammatory response, increasing the expression of other secreted factors like IL-6 and IL-8. IL-1 α was found to affect the cell cycle of osteosarcoma cell to reduce cell growth.


References

1. Dinarello CA. et al. (2009). Annu Rev Immunol. 27:519-50.
2. Dinarello CA. et al. (1998). Int Rev Immunol. 16:457-99.
3. Infante J. et al. (2007). Dement Geriatr Cogn Disord. 23:215-8.
4. Apte RN. et al. (2008). Immunol Rev. 222:222-41.
5. Palmer G. et al. (2005). Cytokine. 2:153-60.
6. Rider P. et al. (2013). Semin Immunol. 25(6):430-8.
7. Orjalo AV. et al. (2009). Proc Natl Acad Sci U S A. 106(40):17031-6.


Synonyms

IL1a, IL1-alpha, Il 1a, IL1 alpha, IL-1, IL1

Selected Validation Data



Purity of Mouse IL-1 alpha was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.



The purity of Mouse IL-1 alpha was greater than 90% as determined by SEC-HPLC.