

# Recombinant Mouse Transferrin protein (His Tag)

Catalog Number: Eg1066

## Basic Information

**Species:**  
Mouse**Purity:**  
>90 %, SDS-PAGE**Tag:**  
His Tag

## Technical Specifications

**Purity:**

&gt;90 %, SDS-PAGE

**Endotoxin Level:**

&lt;0.1 EU/ µg protein, LAL method

**Source:**

HEK293-derived Mouse Transferrin protein Val20-His697 (Accession# Q921I1) with a His tag at the C-terminus.

**GeneID:**

22041

**Accession:**

Q921I1

**Predicted Molecular Mass:**

76.0 kDa

**SDS-PAGE:**

65-80 kDa, reducing (R) conditions

**Formulation:**

Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

## Biological Activity

Not tested

## Storage and Shipping

**Storage:**

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

**Shipping:**

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

## Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

## Background

Transferrin, also known as Trf, is a crucial iron-binding transport protein. It plays a vital role in iron homeostasis by binding and transporting iron throughout the body, particularly from the sites of absorption to the sites of storage and utilization. Transferrin is primarily synthesized by the liver. It can bind two ferric iron (Fe<sup>3+</sup>) ions, each associated with a bicarbonate anion, which helps maintain iron in a soluble form and minimizes its redox activity. Transferrin also interacts with clotting factors and plays a central role in coagulation balance. It has been shown to potentiate the enzymatic activity of thrombin and FXIIa and inhibit antithrombin (AT) independently of iron. Elevated levels of transferrin have been associated with atherosclerosis, and it has been found in increased amounts in the plasma of mice fed a high-fat diet, which correlates with atherosclerotic plaque development.


## References

1. Julia T Bu. et al. (2015). Biometals. 28(3):473-480.
2. Xiaopeng Tang. et al.(2020). Cell Res. 30(2):119-132.
3. Konstantinos Gkouvatsos. et al. (2012). Biochim Biophys Acta. 1820(3):188-202.

## Synonyms

Beta-1 metal-binding globulin, Serotransferrin, Siderophilin, Tf, Trf

## Selected Validation Data



Purity of Recombinant Mouse Transferrin was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.

For technical support and original validation data for this product please contact

T: 027-87531629

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.