

For Research Use Only

Recombinant Mouse CD21 protein (rFc Tag)



Catalog Number: Eg1396

Basic Information	Species: Mouse	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <1.0 EU/ µg protein, LAL method Source: HEK293-derived Mouse CD21 protein Ile12-Trp963 (Accession# P19070) with a rabbit IgG Fc tag at the C-terminus. GeneID: 12902 Accession: P19070 Predicted Molecular Mass: 130.9 kDa SDS-PAGE: Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background			
References			
Synonyms			

Selected Validation Data

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.