

For Research Use Only

Recombinant Human VEGF-C protein (rFc Tag)



Catalog Number: Eg1482

Basic Information

Species:
Human

Purity:
>90 %, SDS-PAGE

Tag:
rFc Tag

Technical Specifications

Purity:

>90 %, SDS-PAGE

Endotoxin Level:

<0.1 EU/ μ g protein, LAL method

Source:

HEK293-derived Human VEGF-C protein Ala112-Arg227 (Accession# P49767) with a rabbit IgG Fc tag at the C-terminus.

GeneID:

7424

Accession:

P49767

Predicted Molecular Mass:

39.4kDa

SDS-PAGE:

42-50 kDa, reducing (R) conditions

Formulation:

Lyophilized from 0.22 μ m filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Not tested

Storage and Shipping

Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

Vascular endothelial growth factor C (VEGF-C) is a protein that is a member of the platelet-derived growth factor / vascular endothelial growth factor (PDGF/VEGF) family. The main function of VEGF-C is to promote the growth of lymphatic vessels (lymphangiogenesis). It acts on lymphatic endothelial cells (LECs) primarily via its receptor VEGFR-3 promoting survival, growth and migration. Apart from vascular targets, VEGF-C is also important for neural development and blood pressure regulation.


References

1. V Joukov. et al. (1996) EMBO J. 1996 Jan 15;15(2):290-98.
2. Barbara Le Bras. et al. (2006) Nat Neurosci. Mar;9(3):340-8.
3. Agnes Machnik. et al. (2009) Nat Med. 15(5):545-52.

Synonyms

VEGFC, vascular endothelial growth factor C, Vascular endothelial growth factor-related protein, VEGF C

Selected Validation Data



Purity of Recombinant Human VEGF-C was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.