

For Research Use Only


Recombinant Mouse VEGF-C protein (rFc Tag)



Catalog Number: Eg1640

Basic Information	Species: Mouse	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/ µg protein, LAL method Source: HEK293-derived Mouse VEGF-C protein Ala108-Arg223 (Accession# P97953-1) with a rabbit IgG Fc tag at the C-terminus. GenelD: 22341 Accession: P97953-1 Predicted Molecular Mass: 39.3 kDa SDS-PAGE: 40-48 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	Vascular endothelial growth factor C (VEGF-C) is a protein that is a member of the platelet-derived growth factor / vascular endothelial growth factor (PDGF/VEGF) family. The main function of VEGF-C is to promote the growth of lymphatic vessels (lymphangiogenesis). It acts on lymphatic endothelial cells (LECs) primarily via its receptor VEGFR-3 promoting survival, growth and migration. Apart from vascular targets, VEGF-C is also important for neural development and blood pressure regulation.		
References	1. V Joukov. et al. (1996) EMBO J. 15(2):290-98. 2. Barbara Le Bras. et al. (2006) Nat Neurosci. 9(3):340-8. 3. Agnes Machnik. et al. (2009) Nat Med. 15(5):545-52.		
Synonyms	Flt4 ligand, Flt4-L, Vascular endothelial growth factor C, Vascular endothelial growth factor-related protein, Vegfc		

Selected Validation Data



Purity of Recombinant Mouse VEGF-C was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) conditions and stained using Coomassie blue.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.