


Recombinant Human B7-H6 protein (rFc Tag) (HPLC verified)


Catalog Number: Eg2568

Basic Information	Species: Human	Purity: >90 %, SDS-PAGE >90 %, SEC-HPLC	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE >90 %, SEC-HPLC Endotoxin Level: <0.1 EU/ µg protein, LAL method Source: HEK293-derived Human B7-H6 protein Asp25-Ser262 (Accession# Q68D85) with a rabbit IgG Fc tag at the C-terminus. GeneID: 374383 Accession: Q68D85 Predicted Molecular Mass: 53.0 kDa SDS-PAGE: 50-70 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	B7-H6, also known as NCR3LG1, is a member of the B7 family of immunoreceptors. It is a 454-aa-long type I transmembrane protein consisting of an extracellular region, a transmembrane segment, and a cytoplasmic domain which contains an immunoreceptor tyrosine-based inhibitory motif (ITIM), an SH2-binding domain, and an SH3-binding motif. B7-H6 triggers Nkp30-mediated activation of human NK cells. It has been reported that B7H6 is expressed in various types of tumor cells and plays a prominent role in natural killer cell immunosurveillance. The apparent molecular weight of B7-H6 is larger than the calculated molecular weight of 51 kDa due to glycosylation.		
References	1. Cameron S Brandt. et al. (2009). J Exp Med. 206(7):1495-503. 2. Alaleh Mohammadi. et al. (2022). Life Sci.304:120709. 3. Ying Chen. et al. (2018). Pathol Oncol Res.4(4):717-721.		
Synonyms	B7H6, NCR3LG1, B7 H6, B7 homolog 6, Natural cytotoxicity triggering receptor 3 ligand 1		

Selected Validation Data



Purity of Recombinant Human B7-H6 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.



The purity of Human B7-H6 was greater than 90% as determined by SEC-HPLC.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.