

For Research Use Only

# Recombinant Mouse Betacellulin/BTC protein (rFc Tag) (HPLC verified)



Catalog Number: Eg2895

## Basic Information

**Species:**  
Mouse

**Purity:**  
>90 %, SDS-PAGE<br>>90 %, SEC-HPLC

**Tag:**  
rFc Tag

## Technical Specifications

**Purity:**  
>90 %, SDS-PAGE<br>>90 %, SEC-HPLC

**Endotoxin Level:**  
<0.1 EU/ μg protein, LAL method

**Source:**  
HEK293-derived Mouse Betacellulin protein Asp32-Gln118 (Accession# Q543J8) with a rabbit IgG Fc tag at the C-terminus.

**GeneID:**  
12223

**Accession:**  
Q543J8

**Predicted Molecular Mass:**  
35.9 kDa

**SDS-PAGE:**  
40-55 kDa, reducing (R) conditions

**Formulation:**  
Lyophilized from 0.22 μm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

## Biological Activity

Not tested

## Storage and Shipping

### Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

### Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

## Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

## Background

Betacellulin (BTC), a member of the epidermal growth factor (EGF) family, binds and activates ErbB1 and ErbB4 homodimers. BTC is synthesized primarily as a transmembrane precursor, which is then processed to mature molecule by proteolytic events. This protein is a ligand for the EGF receptor. BTC was expressed in tumors and involved in tumor growth progression.


## References

- 1.Harris RC, Chung E, Coffey RJ.(2003)Exp Cell Res. 284 (1):2-13.
- 2.Dunbar AJ,Goddard C.(2000) Int J Biochem Cell Biol. 32 (8):805-15.


## Synonyms

Bcn, Betacellulin, Btc, Probetacellulin

## Selected Validation Data



Purity of Recombinant Mouse Betacellulin was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.



The purity of Mouse Betacellulin was greater than 90% as determined by SEC-HPLC.

For technical support and original validation data for this product please contact

T: 027-87531629

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.