

For Research Use Only

Recombinant Mouse ACE protein (rFc Tag)



Catalog Number: Eg3036

Basic Information

Species:
Mouse

Purity:
>90 %, SDS-PAGE

Tag:
rFc Tag

Technical Specifications

Purity:

>90 %, SDS-PAGE

Endotoxin Level:

<0.1 EU/ μ g protein, LAL method

Source:

HEK293-derived Mouse ACE protein Leu35-Gln1264 (Accession# P09470-1) with a rabbit IgG Fc tag at the C-terminus.

GeneID:

11421

Accession:

P09470-1

Predicted Molecular Mass:

167.9 kDa

SDS-PAGE:

140-200 kDa, reducing (R) conditions

Formulation:

Lyophilized from 0.22 μ m filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Not tested

Storage and Shipping

Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background


Angiotensin I Converting Enzyme (ACE) is a transmembrane zinc metallopeptidase that hydrolyzes peptide chains by removing the dipeptide from the C-terminus. ACE is a key component of the renin-angiotensin system (RAS), which plays an important role in blood pressure homeostasis through the production of the vasoconstrictor peptide angiotensin II. ACE also plays a role in fertility, immunity, hematopoiesis, and diseases such as obesity, fibrosis, and Alzheimer's disease dementia. ACE is relevant to immune regulation in that endogenous ACE has been detected in a variety of immune cell populations including macrophages and neutrophils, and its overexpression results in enhanced bactericidal and antitumor responses independent of angiotensin II.

References

1. Oosthuizen D, et al. (2022). J Renin Angiotensin Aldosterone Syst. 2022:9028969.
2. Rahimi Z. (2012). J Nephropathol. 1(3):143-51.
3. Lubbe L, et al. (2020). Clin Sci (Lond). 134(21):2851-2871.

Synonyms

Selected Validation Data



Purity of Recombinant Mouse ACE was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) conditions and stained using Coomassie blue.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.