

For Research Use Only


Recombinant Mouse uPAR/CD87 protein (rFc Tag)



Catalog Number: Eg3156

Basic Information	Species: Mouse	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	<p>Purity: >90 %, SDS-PAGE</p> <p>Endotoxin Level: <0.1 EU/ µg protein, LAL method</p> <p>Source: HEK293-derived Mouse uPAR protein Leu24-Thr297 (Accession# P35456-1) with a rabbit IgG Fc tag at the C-terminus.</p> <p>GeneID: 18793</p> <p>Accession: P35456-1</p> <p>Predicted Molecular Mass: 56.0 kDa</p> <p>SDS-PAGE: 65-90 kDa, reducing (R) conditions</p> <p>Formulation: Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.</p>		
Biological Activity	Not tested		
Storage and Shipping	<p>Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.</p> <ul style="list-style-type: none">• Until expiry date, -20°C to -80°C as lyophilized proteins.• 3 months, -20°C to -80°C under sterile conditions after reconstitution. <p>Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.</p>		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	<p>uPAR/CD87 is a highly glycosylated GPI-anchored membrane protein. In addition to the membrane-anchored form, uPAR is released from the plasma membrane by cleavage of the GPI anchor and can be found as a soluble form (suPAR). uPAR contains three homologous domains (D1-D3) of which the N-terminal one (D1) represents the uPA-binding domain. After binding to uPAR, uPA cleaves plasminogen, generating the active protease plasmin which is involved in a wide variety of physiologic and pathologic processes. In addition to regulating proteolysis, uPAR has important function in cell adhesion, migration and proliferation. Studies reveal that uPAR expression is elevated during inflammation and tissue remodelling and in many human cancers, in which it frequently indicates poor prognosis. suPAR has been detected in plasma, and increased plasma concentrations of suPAR have been found in patients with some advanced cancers.</p>		
References	<ol style="list-style-type: none">1. Ploug M. et al. (1991). J Biol Chem. 266(3):1926-1933.2. Stephens RW. et al. (1997). Clin Chem. 43(10):1868-1876.3. Blasi F. et al. (2002). Nat Rev Mol Cell Biol. 3(12):932-943.4. Smith HW. et al. (2010). Nat Rev Mol Cell Biol. 11(1):23-36.		
Synonyms	CD87, Plaur, uPAR, U-PAR, Urokinase plasminogen activator surface receptor		

Selected Validation Data



Purity of Recombinant Mouse uPAR was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) conditions and stained using Coomassie blue.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.