

For Research Use Only


Recombinant Mouse ATP1B1 protein (mFc Tag)



Catalog Number: Eg3915

Basic Information	Species: Mouse	Purity: >90 %, SDS-PAGE	Tag: mFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/ µg protein, LAL method Source: HEK293-derived Mouse ATP1B1 protein Glu63-Ser304 (Accession# P14094) with a mouse IgG Fc tag at the N-terminus. GeneID: 11931 Accession: P14094 Predicted Molecular Mass: 55.8 kDa SDS-PAGE: 58-70 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	ATP1B1 is one of beta subunits of the Na⁺/K⁺ ATPase and responsible for formation and structural integrity of the Na⁺/K⁺ ATPase. The Na⁺/K⁺ ATPase is a plasma membrane pump consisting of alpha, beta, and gamma subunits. At least four of Na⁺/K⁺-ATPase beta subunits (β 1, β 2, β 3, β 4) have been identified in mammalian cells; the β 1-subunit (ATP1B1) is the most ubiquitous. The Na⁺/K⁺ ATPase β subunits have multiple N-glycosylation sites. The predicted MW of ATP1B1 is 35 kDa, while it migrates around 40-52 kDa due to the variable glycosylation.		
References	1. Cutler CP, et al. (2000). Am J Physiol Regul Integr Comp Physiol. 279(1): R222-9. 2. Liu CG, et al. (2007). Clin Exp Pharmacol Physiol. 34(10):998-1004. 3. Wang J, et al. (2018). Exp Neurol. 302:214-229.		
Synonyms	Atp4b, Sodium/potassium-dependent ATPase subunit beta-1, Sodium/potassium-transporting ATPase subunit beta-1		

Selected Validation Data



Purity of Recombinant Mouse ATP1B1 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.