

Recombinant Human Ephrin B1 protein (mFc Tag) (HPLC verified)

Catalog Number: Eg4188

Basic Information

Species:
Human**Purity:**
>90 %, SDS-PAGE
>90 %, SEC-HPLC**Tag:**
mFc Tag

Technical Specifications

Purity:
>90 %, SDS-PAGE
>90 %, SEC-HPLC**Endotoxin Level:**
<0.1 EU/ µg protein, LAL method**Source:**
HEK293-derived Human Ephrin B1 protein Leu28-Lys237 (Accession# P98172) with a mouse IgG Fc tag at the C-terminus.**GeneID:**
1947**Accession:**
P98172**Predicted Molecular Mass:**
49.6 kDa**SDS-PAGE:**
52-65 kDa, reducing (R) conditions**Formulation:**
Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Not tested

Storage and Shipping

Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

EFNB1 (Ephrin B1) is a member of the ephrin family, which are membrane-bound guidance cues that interact with Eph receptors, a family of receptor tyrosine kinases. EFNB1 is X-inactivated, and it has been proposed that in heterozygous females, patchy loss of ephrin B1 disturbs tissue boundary formation at the developing coronal suture. EFNB1 mainly plays a central part in cell adhesion, angiogenesis, and development of the nervous system.


References

1. Shi Y, et al. (2021). J Oncol. vol. 4701680. Nov.16.
2. Wallis D, et al. (2008). Am J Med Genet A.146A(15):2008-2012.
3. Arthur A, et al. (2019). Exp Hematol. 69:43-53.


Synonyms

EFL3, EFNB1, EPLG2, LERK2, CFND

Selected Validation Data



Purity of Recombinant Human Ephrin B1 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.



The purity of Human Ephrin B1 was greater than 90% as determined by SEC-HPLC.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.