


Recombinant Canine CXCL8/IL-8 protein (rFc Tag)(HPLC verified)


Catalog Number: Eg4552

Basic Information	Species: Canine	Purity: >90 %, SDS-PAGE >90%, SEC-HPLC	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE >90%, SEC-HPLC Endotoxin Level: <0.1 EU/ µg protein, LAL method Source: HEK293-derived Canine CXCL8 protein Ala23-Pro101 (Accession# P41324) with a rabbit IgG Fc tag at the C-terminus. GeneID: 403850 Accession: P41324 Predicted Molecular Mass: 35.4 kDa SDS-PAGE: 36-42 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	Interleukin-8 (IL-8), also known as CXCL8, is a member of the CXC chemokine family. This chemokine is secreted by a variety of cell types including monocytes/macrophages, T cells, neutrophils, fibroblasts, endothelial cells and various tumour cell lines in response to inflammatory stimuli. IL-8 has two main functions. It induces chemotaxis in target cells, primarily neutrophils but also other granulocytes, causing them to migrate to the site of infection. IL-8 also induces phagocytosis once they have arrived. This gene is thought to play a role in the pathogenesis of bronchiolitis, a common respiratory disease caused by viral infection. IL-8 is also known to be a potent promoter of angiogenesis. IL-8 has been implicated in tumour angiogenesis, metastasis and poor prognosis in breast cancer. IL-8 may represent a novel therapeutic target for estrogen-driven breast carcinogenesis and tumour progression.		
References	1.Ha H. et al. (2017) Theranostics. 7(6):1543-1588. 2.Matsushima K. et al. (2022) Cytokine. 153:155828. 3.Todorović-Raković N. et al. (2013) J Interferon Cytokine Res. 33(10):563-70. 4.Tatsuno R. et al. (2024) Cell Death Dis. 15(7):471.		
Synonyms	CXCL8, IL-8, IL8, Interleukin-8		

Selected Validation Data



Purity of Recombinant Canine CXCL8 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.



The purity of Canine CXCL8/IL-8 was greater than 90% as determined by SEC-HPLC.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.