

For Research Use Only


Recombinant Human IFNA10 protein (rFc Tag)



Catalog Number: Eg4952

Basic Information	Species: Human	Purity: >90 %, SDS-PAGE	Tag: rFc Tag
Technical Specifications	Purity: >90 %, SDS-PAGE Endotoxin Level: <0.1 EU/ µg protein, LAL method Source: HEK293-derived Human IFNA10 protein Cys24-Asp189 (Accession# P01566) with a rabbit IgG Fc tag at the N-terminus. GeneID: 3446 Accession: P01566 Predicted Molecular Mass: 46.6 kDa SDS-PAGE: 45-52 kDa, reducing (R) conditions Formulation: Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.		
Biological Activity	Not tested		
Storage and Shipping	Storage: It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. • Until expiry date, -20°C to -80°C as lyophilized proteins. • 3 months, -20°C to -80°C under sterile conditions after reconstitution. Shipping: The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.		
Reconstitution	Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.		
Background	Interferons (IFN) are a family of cytokines with potent anti viral, antiproliferative and immunomodulatory properties, classified based on their binding specificity to cell surface receptors. IFNA10 has been shown to be a strong inducer of IFN-stimulated genes and anti viral protection.		
References	1. Pestka S. et al. (1987). Annu. Rev. Biochem. 56:727. 2. Moll H.P. et al. (2011). Cytokine. 53:52.		
Synonyms	IFN-alpha-10, Interferon alpha-10, Interferon alpha-6L, Interferon alpha-C, LeIF C		

Selected Validation Data



Purity of Recombinant Human IFNA10 was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.