

Recombinant Human Podoplanin protein (rFc Tag)(HPLC verified)

Catalog Number: Eg5024

Basic Information

Species:
Human**Purity:**
>90 %, SDS-PAGE
>90 %, SEC-HPLC**Tag:**
rFc Tag

Technical Specifications

Purity:
>90 %, SDS-PAGE
>90 %, SEC-HPLC**Endotoxin Level:**
<0.1 EU/ µg protein, LAL method**Source:**
HEK293-derived Human Podoplanin protein Ala23-Leu131 (Accession# Q86YL7-1) with a rabbit IgG Fc tag at the C-terminus.**GeneID:**
10630**Accession:**
Q86YL7-1**Predicted Molecular Mass:**
37.4 kDa**SDS-PAGE:**
48-68 kDa, reducing (R) conditions**Formulation:**
Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Not tested

Storage and Shipping

Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

Podoplanin (PDPN) is a transmembrane glycoprotein that plays a role in a variety of biological processes including regulation of blood-lymphatic vessel development, cerebral vascular pattern and integrity, cell motility, tumorigenesis and metastasis, and mediates tumor cell-induced platelet aggregation in different cancer types. It is a lymphatic marker because the expression of podoplanin has been detected in lymphatic but not blood vascular endothelium, and is useful as the marker of tumor-associated Lymphangiogenesis. When expressed in keratinocytes, PDPN induces changes in cell morphology with transfected cells showing an elongated shape, numerous membrane protrusions, major reorganization of the actin cytoskeleton, increased motility and decreased cell adhesion.


References

1. Ukaji T, et al. Cancer Sci. 2021; 112(6):2299-2313.
2. Krishnan H, et al. J Biol Chem. 2013; 288(17):12215-21.
3. Wang X, et al. Cancer Gene Ther. 2023; 30(2):345-357.
4. Astarita JL, et al. Front Immunol. 2012; 3:283.


Synonyms

PDPN, 29kDa cytosolic podoplanin intracellular domain, Aggrus, PA2.26 antigen, PSEC0003

Selected Validation Data



Purity of Recombinant Human Podoplanin was determined by SDS-PAGE. The protein was resolved in an SDS-PAGE in reducing (R) and non-reducing (NR) conditions and stained using Coomassie blue.



The purity of Human Podoplanin was greater than 90% as determined by SEC-HPLC.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.