

Catalog Number:

Basic Information

Technical Specifications

Purity:
Endotoxin Level:
Source:
GeneID:
Accession:
Predicted Molecular Mass:
SDS-PAGE:
Formulation:

Biological Activity

Storage and Shipping

Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Background

References

Synonyms

Selected Validation Data

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgc.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.