

For Research Use Only

Recombinant Human SLITRK6 protein (rFc Tag)



Catalog Number: Eg5865

Basic Information

Species:
Human

Purity:
>90 %, SDS-PAGE

Tag:
rFc Tag

Technical Specifications

Purity:

>90 %, SDS-PAGE

Endotoxin Level:

<0.1 EU/ μ g protein, LAL method

Source:

CHO-derived Human SLITRK6 protein Ser27-Ser608 (Accession# Q9H5Y7) with a rabbit IgG Fc tag at the C-terminus.

GeneID:

84189

Accession:

Q9H5Y7

Predicted Molecular Mass:

90.4 kDa

SDS-PAGE:

Formulation:

Lyophilized from 0.22 μ m filtered solution in 0.5M arg, 30mM Tris, pH8.0. Normally 5% trehalose and 5% mannitol are added as protectants before lyophilization.

Biological Activity

Not tested

Storage and Shipping

Storage:

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

- Until expiry date, -20°C to -80°C as lyophilized proteins.
- 3 months, -20°C to -80°C under sterile conditions after reconstitution.

Shipping:

The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature.

Reconstitution

Briefly centrifuge the tube before opening. Reconstitute at 0.1-0.5 mg/mL in sterile water.

Background

References

Synonyms

Selected Validation Data

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.