

For Research Use Only

# FITC Anti-Mouse CSF1R/CD115 (AFS98)

Catalog Number: FITC-65095



## Basic Information

Catalog Number:

FITC-65095

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC036343

GeneID (NCBI):

12978

UNIPROT ID:

P09581

Full Name:

colony stimulating factor 1 receptor

Purification Method:

Affinity purification

CloneNo.:

AFS98

Excitation/Emission maxima  
wavelengths:

498 nm / 526 nm

## Applications

Tested Applications:

FC

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

Phosphate based buffer with 0.09% sodium azide and 0.1% gelatin, pH 7.2.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> Balb/c mouse peritoneal macrophages were surface stained with 0.5 ug FITC Anti-Mouse CD115 (c-fms) (FITC-65095, Clone: AFS98) (green) or 0.5 ug FITC-Rat IgG2a Isotype Control (black). Cells were not fixed.