

For Research Use Only

FITC Plus Anti-Mouse CD11c (N418)

Catalog Number: FITC-65130

6 Publications



Basic Information

Catalog Number:

FITC-65130

Size:

100ug, 0.5 mg/mL

Source:

Armenian Hamster

Isotype:

IgG

GenBank Accession Number:

BC167225

GeneID (NCBI):

16411

UNIPROT ID:

Q9QXH4

Full Name:

integrin alpha X

Purification Method:

Affinity purification

CloneNo.:

N418

Excitation/Emission maxima wavelengths:

495 nm / 524 nm

Applications

Tested Applications:

FC

Cited Applications:

IF, FC

Species Specificity:

Mouse

Cited Species:

mouse

Background Information

Integrins are cell adhesion receptors that are heterodimers composed of non-covalently associated α and β subunits (PMID: 9779984). CD11c, also known as integrin αX , is a type I transmembrane glycoprotein present on a variety of cells, including monocytes/macrophages, granulocytes, a subset of B cells, NK cells and dendritic cells (PMID: 2897326; 1680915; 1694698; 17389580). As a result of its high level of expression on most dendritic cells, CD11c is typically considered to be a marker of conventional dendritic cells (PMID: 27119555). CD11c forms an α / β heterodimer with CD18 (integrin $\beta 2$). CD11c/CD18 acts a receptor for fibrinogen and is important in monocyte adhesion and chemotaxis (PMID: 1671533).

Notable Publications

Author	Pubmed ID	Journal	Application
Lian Zhu	36321642	J Mater Chem B	FC
Ze Wang	38883315	Bioact Mater	FC
Yi Chen	38768811	J Adv Res	IF

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface stained with APC Anti-Mouse MHC Class II (I-A/I-E) (APC-65122, Clone: M5/114.15.2) and 0.25 ug FITC Plus Anti-Mouse CD11c (FITC-65130, Clone: N418) or 0.25 ug FITC Plus Hamster IgG isotype control. Cells were not fixed.