

For Research Use Only

# FITC Plus Anti-Human STRO-1 (STRO-1)



Catalog Number: FITC-65184

1 Publications

## Basic Information

Catalog Number:

FITC-65184

Size:

100 tests, 5 µl/test

Source:

Mouse

Isotype:

IgM, lambda

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Purified by precipitation and chromatography

CloneNo.:

STRO-1

Excitation/Emission maxima wavelengths:

495 nm / 524 nm

## Applications

Tested Applications:

FC

Cited Applications:

FC

Species Specificity:

Human

Cited Species:

human

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Tao Jiang	37441553	J Tissue Eng	FC

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> human whole blood cells were surface stained with 5 ul FITC Plus Anti-Human STRO-1 (FITC-65184, Clone: STRO-1) or unstained. Cells were not fixed. CD45dim, CD34+ cells were gated.