

For Research Use Only

FITC Plus Anti-Rat CD45 (OX1)

Catalog Number: FITC-65228

1 Publications



Basic Information

Catalog Number:

FITC-65228

Size:

100ug, 500 µg/mL

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

XM_006249910

GeneID (NCBI):

24699

UNIPROT ID:

P04157

Full Name:

protein tyrosine phosphatase,
receptor type, C

Calculated MW:

143 kDa

Purification Method:

Affinity purification

CloneNo.:

OX1

Excitation/Emission maxima
wavelengths:

495 nm / 524 nm

Applications

Tested Applications:

FC

Cited Applications:

FC

Species Specificity:

rat

Cited Species:

rat

Background Information

CD45, also known as protein tyrosine phosphatase, receptor type C, is a type I transmembrane protein expressed on the surface of all haematopoietic cells with the exception of erythrocytes and platelets (PMID: 3489673; 28615666). CD45 is a pan-haematopoietic cell marker and has been shown to be essential for T- and B-cell activation and signalling (PMID: 9429890; 16378097).

Notable Publications

Author	Pubmed ID	Journal	Application
Yu Meng	39339434	Molecules	FC

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ Wistar rat splenocytes were surface stained with 1 ug FITC Plus Anti-Rat CD45 (FITC-65228, Clone: OX1) or 1 ug Control Antibody. Cells were not fixed.