

For Research Use Only

FITC Plus Anti-Mouse CD94 (18D3)

Catalog Number: FITC-65240



Basic Information

Catalog Number:

FITC-65240

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgG2a, kappa

GenBank Accession Number:

BC120854

GeneID (NCBI):

16643

UNIPROT ID:

O54707

Full Name:

killer cell lectin-like receptor,
subfamily D, member 1

Purification Method:

Affinity purification

CloneNo.:

18D3

Excitation/Emission maxima
wavelengths:

495 nm / 524 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface co-stained with CoraLite®647 Anti-Mouse CD49b (HMA2) (CL647-65245, Clone: HMA2) and 0.5 ug FITC Plus Anti-Mouse CD94 (FITC-65240, Clone: 18D3) or 0.5 ug Isotype Control. Cells were not fixed.