

For Research Use Only

# FITC Plus Mouse IgM Isotype Control (MOPC-104E)



Catalog Number: FITC-65243

## Basic Information

Catalog Number:

FITC-65243

Size:

100ug, 200 µg/ml

Source:

Mouse

Isotype:

IgM, lambda

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

N/A

CloneNo.:

MOPC-104E

Excitation/Emission maxima  
wavelengths:

495 nm / 524 nm

## Applications

Tested Applications:

FC

Species Specificity:

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1X10<sup>6</sup> human whole blood cells were surface stained with 0.2 ug FITC Plus Anti-Mouse IgM Isotype Control (red), 5 ul FITC Plus Anti-Human STRO-1 (green), or unstained (blue). Cells were not fixed. CD45dim, CD34 lymphocytes were gated.