

For Research Use Only

# FcZero-rAb® FITC Plus Rabbit IgG Isotype Control Recombinant Antibody

Catalog Number: FITC-FcA98136



## Basic Information

Catalog Number:

FITC-FcA98136

Concentration:

100ug, 500 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

240953C9

Recommended Dilutions:

FC: 0.25 ug per 10<sup>6</sup> cells in a 100 µl suspension

Excitation/Emission maxima wavelengths:

498 nm / 526 nm

## Applications

Tested Applications:

FC

Species Specificity:

n/a

Positive Controls:

FC : human PBMCs,

## Background Information

This product is an isotype control antibody, which is used to estimate the non-specific binding of target primary antibodies due to Fc receptor binding or other protein-protein interactions. An isotype control antibody should have the same immunoglobulin type and be used at the same concentration as the test antibody.

## Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1x10<sup>6</sup> human PBMCs were unstained (left, blue in histogram), or surface stained with 0.25 ug FITC Rabbit IgG Isotype Control RecAb (FITC-FcA98136, Clone: 240953C9) (mid, red in histogram) or 0.25 ug FITC Plus Anti-Human CD14 Rabbit RecAb (FITC-FcA98040, Clone: 230332D7) (right, green in histogram). Cells were not fixed.